



**Venous Ulcer Compression System AVAILABLE ON DRUG TARIFF FP10/GP10**

**It is essential that all sections are completed fully; otherwise your order cannot be processed**

**PHARMACY/HOSPITAL TO COMPLETE**


Order Date: \_\_\_\_\_ Order No: \_\_\_\_\_ Invoice to: Wholesaler  Buying Group   
 Name and Address of Pharmacy/Hospital: \_\_\_\_\_ Name of Wholesaler/Buying Group: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_ Postcode: \_\_\_\_\_ Account Number: \_\_\_\_\_  
 Telephone No: \_\_\_\_\_ Patient's Name: \_\_\_\_\_

**Measuring Guide**

	Small	Medium	Large	X Large	XX Large
B: Calf	32 - 38	34 - 40	36 - 42	38 - 45	40.5 - 50
C: Ankle	19 - 25.5	21 - 27.5	23 - 29.5	25 - 31.5	26.5 - 34

Measurements in cm

Measurements should be taken with the person seated, and feet flat on the floor. If the person has a swollen leg, measurements should be taken first thing in the morning. <sup>1</sup>



Leg Length (L)	
Short	34 - 38
Regular	38 - 42
Long	42 - 46

References: <sup>1</sup> NICE Clinical Knowledge Summaries



Altipress 40<sup>®</sup> Hosiery Kit (1x Class III Soft Beige Stocking and 2 x Liners)

Kit	Length	Small (PIP Code)	Qty	Medium (PIP Code)	Qty	Large (PIP Code)	Qty	X Large (PIP Code)	Qty	XX Large (PIP Code)	Qty
Altipress 40 <sup>®</sup> Hosiery Kit	Short	337-1325		337-1333		337-1341		337-1358		337-1366	
Altipress 40 <sup>®</sup> Hosiery Kit	Regular	337-1374		337-1382		337-1390		337-1408		337-1416	
Altipress 40 <sup>®</sup> Hosiery Kit	Long	337-1424		337-1432		337-1440		337-1457		337-1465	



Altipress<sup>®</sup> Liner Pack (3 Liners)

Kit	Length	Small (PIP Code)	Qty	Medium (PIP Code)	Qty	Large (PIP Code)	Qty	X Large (PIP Code)	Qty	XX Large (PIP Code)	Qty
<b>White</b>											
Altipress <sup>®</sup> Liner Pack	Short	337-1473		337-1481		337-1499		337-1507		337-1515	
Altipress <sup>®</sup> Liner Pack	Regular	337-1523		337-1531		337-1549		337-1556		337-1564	
Altipress <sup>®</sup> Liner Pack	Long	337-1572		337-1580		337-1598		337-1606		337-1614	
<b>Soft Beige</b>											
Altipress <sup>®</sup> Liner Pack	Short	341-8506		341-8696		341-8688		341-8670		341-8613	
Altipress <sup>®</sup> Liner Pack	Regular	341-8662		341-8654		341-8647		341-8621		341-8605	
Altipress <sup>®</sup> Liner Pack	Long	341-8598		341-8563		341-8530		341-8522		341-5414	

**Fax your order to 01509 501721 or email to enquiries@altimed.co.uk**

For additional forms visit [www.altimed.co.uk](http://www.altimed.co.uk), email [enquiries@altimed.co.uk](mailto:enquiries@altimed.co.uk) or phone 01509 501720

AltiMed Ltd, an Urgo Medical company. Sullington Road, Shephed, Loughborough LE12 9JG

