

**72 HOUR
DESPATCH***

British Standard AVAILABLE ON DRUG TARIFF FP10/GP10

It is essential that all sections are completed fully; otherwise your order cannot be processed

PHARMACY/HOSPITAL TO COMPLETE

Order Date: _____ Order No: _____ Invoice to: Hospital/Clinic Buying Group
 Wholesaler Pharmacy Chain
 Name and Address of Pharmacy/Hospital: _____ Name of Wholesaler/Buying Group/Pharmacy Chain: _____
 _____ Acc. No. _____
 _____ Postcode: _____ Patient's Name: _____
 Telephone No: _____ Patient's Previous Reference No. (if known): _____

Measurements completed by:

Name: _____ Position: _____
 Healthcentre: _____ Postcode: _____
 Contact Tel No: _____

Compression

Class I (14-17mmHg)
 Class II (18-24mmHg)
 Class III (25-35mmHg)

Style

Below Knee Thigh Length Hold Ups - Plain
 Thigh Length Standard Thigh Length Hold Ups - Lace
Foot: Closed toe Open toe

Quantity (singles)


Left Leg: _____ Right Leg: _____
 Beige Black Beige Black

Measurements should be taken with the person seated, and feet flat on the floor. If the person has a swollen leg, measurements should be taken first thing in the morning¹.

COMMENTS:

Altimed cannot be held responsible for incorrect measurements

Below Knee		
Please fill in all the appropriate measurements in CENTIMETRES	Left	Right
CIRCUMFERENCE CD (under knee)		
CC (calf)		
CB (ankle)		
CA (Measured at base of little toe)		
LENGTH LD (below knee)		
LA (open toe only)		
LZ (closed toe only)		



Thigh Length		
Please fill in all the appropriate measurements in CENTIMETRES	Left	Right
CIRCUMFERENCE CG (thigh)		
CD (under knee)		
CC (calf)		
CB (ankle)		
CA (Measured at base of little toe)		
LENGTH LT (waist to floor)		
LG (leg length)		
LA (open toe only)		
LZ (closed toe only)		

*Altimed Ltd will despatch orders within 72 hours of receipt, subject to forces beyond its control

ALTIFORM 10/13

References: ¹ NICE Clinical Knowledge Summaries.

Fax your order to 01509 501 721 or email to enquiries@altimed.co.uk

For additional forms visit www.altimed.co.uk, email enquiries@altimed.co.uk or phone 01509 501720

Altimed Ltd, an Urgo Medical company. Sullington Road, Shepshed, Loughborough LE12 9JG